



	Grid Pad 8	Grid Pad 10
Operating system	Windows 10 Pro	Windows 10 Enterprise
Software	Grid 3	Grid 3
Processor	Atom Z3795 1.6GHz (2.4GHz Burst)	Atom Z3795 1.6GHz (2.4GHz Burst)
Memory / RAM	4GB	4GB
Hard drive	64GB eMMC	64GB eMMC
Display	8 inch 1920 x 1080	10 inch 1920 x 1200
Weight	1.2Kg	1.4Kg
Battery life	4 hours	7 hours
Charge time	4 hours	5 hours
Tablet camera	Front camera: 2mp Rear camera: 5mp	Front camera: 2mp Rear camera: 5mp
Dimensions	24.5 x 16.5 x 4.5 cm	28.5 x 19.2 x 4.2 cm
Eye gaze	none	none
Access	2 x 3.5mm switch ports 3.5mm Auditory feedback port 3.5mm Headphone port 2 x USB 2.0 port Keyguard compatible	2 x 3.5mm switch ports 3.5mm Auditory feedback port 3.5mm Headphone port 2 x USB 2.0 port Keyguard compatible
Environment control	GEWA Infrared EasyWave or Z-Wave radio	GEWA Infrared EasyWave or Z-Wave radio
Connectivity	Dual band wireless WLAN 802.11 a/b/g/n Bluetooth 4.0 Micro SD Card slot	Dual band wireless WLAN 802.11 a/b/g/n Bluetooth 4.0 Micro SD Card slot
Mounting plate	Daessy and Rehadapt	Daessy and Rehadapt
Other features	Built in stand	Built in stand