



The all-new Grid Pad is designed for communication all day, every day. A long 15-hour battery life, rugged design and multiple access options ensure you always have your voice.

Operating system	Windows 10 Pro
Software	Grid 3
Processor	Intel® Core™ i5-7Y54 Processor, 4MB Cache, up to 3.20 GHz
Memory / RAM	8GB DDR3
Hard drive	SSD SATA III, 128GB
Display	12.5" Hardened Glass Intel HD 1920 X 1080p
Weight	2.7kg 5.9 lbs
Battery life	15 hours typical use 10 hours with eye gaze
Charge time	10 hours
Tablet camera	2MP front, 5MP rear
Dimensions	355 x 223 x 45 mm 14 x 8.8 x 1.8 inch
Eye gaze	Optional Grid Pad camera unit (EyeTech TM5 Mini or Irisbond Duo) Eye gaze mounting plates also available (PCEye Mini, Alea, LC Edge and MyGaze)
Access	2 x 3.5mm switch ports, 3.5mm headphone / mic combo port, touch screen, magnetic Keyguard compatible, 1 x USB 2.0 port, 1 x USB3.0 port (Eye Gaze uses the USB 3.0 port)
Environment control	GEWA Infrared, EasyWave / Z-Wave radio
Connectivity	Dual band Wireless (AC 8265), Bluetooth 4.2 LE
Mounting plate	Daessy and Rehadapt
Warranty	3 years SmartCare (UK only)
Other features	Built in stand and Remote Power Button (batteries included)